

Default Canto Configuration

Laser	Detector Position	Dichroic Filter	Band Pass Filter	Parameter(s)
Blue (488nm)	A	735LP	780/60	PE-Cy7
	B	655LP	670LP	PE-Cy5.5
	C	610LP	-	
	D	556LP	585/42	PE
	E	502LP	530/30	A488, FITC, GFP, YFP, FDG, CFSE
	F	-	488/10	SSC
	G	-	-	
	H	-	-	
Red (633nm)	A	735LP	780/60	A750, APC-Cy7, e780
	B	685LP	-	-
	C	-	660/20	A633, APC, A647, Cy5
Violet (405nm)	A	502LP	510/50	BV510, Ametrine
	B	-	450/50	BV421, DAPI

***Other configurations available. For more information contact the flowcore@scripps.edu**